

<b>Mutant gene</b>	<b>Forward Primer (FP)</b>	<b>Reverse Primer (RP)</b>
<i>COX11-C111A-HA</i>	CGGTACCACTCTATAGAGCC ATTGCTGCTCGTACTGG	CCAGTACGAGCAGCAATGGCTCTA TAGAGTGGTACCG
<i>COX11-C208A-HA</i>	GAAGCAGCTCAATACTTTAA TAAAATCCAAGCTTTTGGCTT TGAAGAACAAAAGCTA	TAGCTTTTGGTCTTCAAAGCAAAA AGCTTGGATTTTATTAAAGTATTG AGCTGCTTC
<i>COX11-C210A-HA</i>	GAAGCAGCTCAATACTTTAA TAAAATCCAATGTTTGCCTT TGAAGAACAAAAGCTAG	CTAGCTTTTGGTCTTCAAAGGCAA AACATTGGATTTTATTAAAGTATT GAGCTGCTTC

**Supplementary Table 2.** The primers used for mutagenesis are shown. Quikchange Lightning mutagenesis kit was used and the pRS416 plasmid containing the corresponding wildtype gene was used as a template for amplification.